

Chill tray Metos SCSL2 with SCBC-2 GN1/2 brass cover

Put the Chill Tray in the freezer overnight (+12 hours). It will maintain a temperature below +8 °C for at least 3 hours..

Material: aluminium with brass cover



Chill tray Metos SCSL2 with SCBC-2 GN1/2 brass cover

| | |
|----------------------------|------------|
| Item width mm | 325 |
| Item depth mm | 265 |
| Item height mm | 36 |
| Package volume | 0.008 |
| Unit of volume | m3 |
| Package volume | 0.008 m3 |
| Package length | 37 |
| Package width | 31 |
| Package height | 7 |
| Package unit of dimension | cm |
| Package dimensions (LxWxH) | 37x31x7 cm |
| Net weight | 2.5 |
| Net weight | 2.5 kg |
| Gross weight | 3 |
| Package weight | 3 kg |
| Unit of weight | kg |